

Anti-Calpain antibody [Calp3C/12A2] (STJ16101689)

STJ16101689

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Calpain is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Rabbit/Mouse/Dog/Chicken/Hamster/Pig/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Calp3C/12A2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:100
Formulation	Lyophilized
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide containing amino acids 355-370 of the human calpain 3 sequence.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081