

Anti-Calpain antibody [Calp3D/2C4] (STJ16101690)

STJ16101690

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Calpain is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Rabbit/Dog/Hamster/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Calp3D/2C4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:25-1:50
Formulation	Lyophilized
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide containing amino acids 1–19 of the human calpain 3 sequence.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081