

Anti-Caveolin-1 antibody [4D6] (STJ16101692)

STJ16101692

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Caveolin-1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4D6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:200 WB
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	857
Gene Symbol	CAV1
Uniprot ID	CAV1_HUMAN
Immunogen	Prokaryotic recombinant protein corresponding to the N-terminus domain of the human caveolin-1 molecule
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081