

Anti-NCAM antibody [CD564] (STJ16101712)

STJ16101712

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD56 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CD564
Concentration	>=11 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:100
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein corresponding to a region of the extracellular domain of the human CD56 molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	

