

Anti-Cytokeratin 20 antibody [K2 20.8] (STJ16101720)

STJ16101720

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cytokeratin 20 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	K2 20.8
Concentration	>=9, 26 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:50
Formulation	PBS (pH 7.6) with 1% BSA and sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Cytoskeletal preparation isolated from microdissected villi of human duodenal mucosa.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081