

Anti-Cytokeratin (5/6/8/18) antibody [5D3 and LP34] (STJ16101721)

STJ16101721

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cytokeratin (5/6/8/18) is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-P/IHC-F
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5D3 and LP34
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:50-1:100 IHC-F 1:50-1:100
Formulation	Tissue culture supernatant with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	5D3, Cytokeratins for the human breast carcinoma cell line MCF-7. LP34, Detergent-insoluble fraction of psoriatic human epidermis.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081