

Anti-Dystrophin antibody (Rod Domain) [DY4/6D3] (STJ16101729)

STJ16101729

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Dystrophin (Rod Domain) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Mouse
Reactivity	Human/Mouse/Rat/Rabbit/Dog/Hamster/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DY4/6D3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Lyophilized
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Bacterial fusion protein
Immunogen Region	Rod Domain
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081