

Anti-Epidermal Growth Factor Receptor antibody [EGFR.113] (STJ16101734)

STJ16101734

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Epidermal Growth Factor Receptor is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EGFR.113
Concentration	>=26 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:20
Formulation	Tissue culture supernatant with Sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1956
Gene Symbol	EGFR
Uniprot ID	EGFR_HUMAN
Immunogen	Prokaryotic fusion protein corresponding to the external domain of the epidermal growth factor receptor molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	