

Anti-Fascin antibody [IM20] (STJ16101740)

STJ16101740

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Fascin is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IM20
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:200-1:400 WB 1:100-1:200
Formulation	Tissue culture supernatant with 15 mM Sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6624
Gene Symbol	FSCN1
Uniprot ID	FSCN1_HUMAN
Immunogen	Prokaryotic recombinant protein corresponding to the C-terminal region of the fascin molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081