

Anti-Gamma-Sarcoglycan antibody [35DAG/21B5] (STJ16101745)

STJ16101745

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Gamma-Sarcoglycan is suitable for use in Immunohistochemistry research applications. |
| Applications | IHC-F |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|------------------------------|
| Clonality | Monoclonal |
| Clone ID | 35DAG/21B5 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | Lyophilized |
| Isotype | IgG2bk |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 6445 |
| Gene Symbol | SGCG |
| Uniprot ID | SGCG_HUMAN |
| Immunogen | Synthetic peptide containing amino acids 167-178 of the rabbit gamma-sarcoglycan sequence. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081