

Anti-IgA antibody [N1CLA] (STJ16101750)

STJ16101750

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Immunoglobulin A is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N1CLA
Concentration	>=46 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:300
Formulation	Tissue culture supernatant with Sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein corresponding to 129 amino acids of the heavy chain constant region of the human immunoglobulin A molecule
Immunogen Region	
Specificity	
Immunogen Sequence	

