

Anti-Interleukin 6 antibody [10C12] (STJ16101753)

STJ16101753

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin 6 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10C12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:50
Formulation	Tissue culture supernatant with 15mM Sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A prokaryotic recombinant protein containing part of the N-terminal domain of the interleukin 6 molecule
Immunogen Region	
Specificity	
Immunogen Sequence	

