

Anti-Kappa Light Chain antibody [CH15] (STJ16101755)

STJ16101755

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Kappa Light Chain is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CH15
Concentration	≥ 125 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:200
Formulation	Tissue culture supernatant with 15mM Sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein corresponding to 106 amino acids of the human kappa light chain molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	

