

Anti-Ki67 Antigen antibody [MM1] (STJ16101756)

STJ16101756

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ki67 antigen is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MM1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tissue culture supernatant with 15mM Sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant fusion protein corresponding to a 1086 bp Ki67 motif-containing cDNA fragment.
Immunogen Region	
Specificity	
Immunogen Sequence	

