

Anti-Lambda Light Chain antibody [SHL53] (STJ16101757)

STJ16101757

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lambda Light Chain is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SHL53
Concentration	>=554 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:200
Formulation	Tissue culture supernatant with 15mM Sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein corresponding to 105 amino acids of the human lambda light chain molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081