

Anti-Merosin Laminin Alpha 2 Chain antibody [Mer3/22B2] (STJ16101762)

STJ16101762

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Merosin Laminin Alpha 2 Chain is suitable for use in Immunohistochemistry research applications.
Applications	IHC-F
Host/Source	Mouse
Reactivity	Human/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Mer3/22B2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Lyophilized
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence