

Anti-p63 antibody [7JUL] (STJ16101780)

STJ16101780

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-p63 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	7JUL
Concentration	>=208 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:25
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant fusion protein corresponding to a region (amino acids 319-410)
Immunogen Region	
Specificity	
Immunogen Sequence	