

Anti-Progesterone Receptor (A/B Forms) antibody [16+SAN27] (STJ16101785)

STJ16101785

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Progesterone Receptor (A/B Forms) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	16+SAN27
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein corresponding to the N-terminal region of the A form of the human progesterone receptor generating clone 16 and a prokaryotic recombinant protein corresponding to the 164 amino acid N-terminal region unique to the B form
Immunogen Region	
Specificity	
Immunogen Sequence	

