

Anti-Thyroid Peroxidase antibody [AC25] (STJ16101796)

STJ16101796

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Thyroid Peroxidase is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AC25
Concentration	>=10.8 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:80
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	7173
Gene Symbol	TPO
Uniprot ID	PERT_HUMAN
Immunogen	Prokaryotic recombinant protein containing 105 amino acids corresponding to a portion of the N-terminus of the human thyroid peroxidase molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081