

Anti-Tyrosinase antibody [T311] (STJ16101798)

STJ16101798

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tyrosinase is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	T311
Concentration	>=89 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:50
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7299
Gene Symbol	TYR
Uniprot ID	TYRO_HUMAN
Immunogen	Recombinant prokaryotic protein corresponding to the tyrosinase molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081