

Anti-CEA-2 antibody [12-140-10] (STJ16101814)

STJ16101814

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CEA-2 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	12-140-10
Concentration	>=38 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:200
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	CEA isolated from liver metastasis of colorectal carcinomas by PCA extraction followed by ion exchange and gel filtration chromatography.
Immunogen Region	
Specificity	
Immunogen Sequence	