

Anti-Dystrophin antibody [13H6] (STJ16101819)

STJ16101819

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Dystrophin is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	13H6
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:10-1:20
Formulation	Lyophilized
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein corresponding to a region of the rod domain of the human dystrophin molecule
Immunogen Region	
Specificity	
Immunogen Sequence	