

Anti-Glial Fibrillary Acidic Protein antibody [GA5] (STJ16101820)

STJ16101820

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Glial Fibrillary Acidic Protein is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GA5
Concentration	>=70 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:400
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2670
Gene Symbol	GFAP
Uniprot ID	GFAP_HUMAN
Immunogen	Porcine spinal cord.
Immunogen Region	
Specificity	
Immunogen Sequence	