

Anti-CM5 antibody (STJ16101825)

STJ16101825

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-p53 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	IHC-P 1:200 WB 1:2000-1:4000
Formulation	PBS with 1% BSA and sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant mouse p53 protein
Immunogen Region	
Specificity	
Immunogen Sequence	