

Anti-Kip2 antibody [25B2] (STJ16101826)

STJ16101826

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-p57 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	25B2
Concentration	>=19 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:50
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant antigen corresponding to a 116 amino acid region of the N-terminus of the p57 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	