

Anti-PTH antibody [105G7] (STJ16101828)

STJ16101828

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Parathyroid Hormone is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	105G7
Concentration	>=12 mg/L
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P 1:300
Formulation	Tissue culture supernatant with sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Prokaryotic recombinant protein corresponding to the entire human parathyroid hormone molecule.
Immunogen Region	
Specificity	
Immunogen Sequence	