

Anti-MDR3 P-glycoprotein antibody [P3II-26] (STJ16101833)

STJ16101833

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MDR3 P-glycoprotein is suitable for use in Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	ICC/IHC-F/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	P3II-26
Concentration	100 ug/mL
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	ICC 1:20-1:50 IHC-F 1:2 WB
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	MDR3 P-gp aa 629-692
Immunogen Region	
Specificity	
Immunogen Sequence	