

Anti-PROSTEIN (P501S) antibody [ZR9] (STJ16101835)

STJ16101835

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-PROSTEIN (P501S) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR9
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:50
Formulation	Tissue culture supernatant with 0.2% BSA and 15mM sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthesized peptides to the N-terminus of human prostein
Immunogen Region	
Specificity	
Immunogen Sequence	