

Anti-Progesterone Receptor antibody [ZR4] (STJ16101837)

STJ16101837

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Progesterone Receptor is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR4
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:200-1:400
Formulation	Tissue culture supernatant with 0.2% BSA and 15mM sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant protein encoding human progesterone receptor 412-526 aa
Immunogen Region	
Specificity	
Immunogen Sequence	