

Anti-Estrogen Receptor antibody [ZR2] (STJ16101842)

STJ16101842

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Estrogen Receptor is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ZR2
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:300
Formulation	Tissue culture supernatant with 0.2% BSA and 15mM sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A synthetic peptide corresponding to residues on the C-terminus of human ER-Alpha
Immunogen Region	
Specificity	
Immunogen Sequence	