

Anti-CD221 IGF Receptor antibody [1H7] (STJ16101843)

STJ16101843

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD221 IGF Receptor is suitable for use in Immunohistochemistry, Flow Cytometry and Western Blot research applications.
Applications	IHC-P/FC/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1H7
Concentration	
Conjugation	Unconjugated
Purification	Purified
Dilution Range	IHC-P 1:10-1:40 FC 1:10 WB
Formulation	PBS with azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified human placental CD221 IGF-1 receptor
Immunogen Region	
Specificity	
Immunogen Sequence	