

## Anti-MART-1 (MELAN A) antibody [M2-9E3] (STJ16101844)

STJ16101844

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-MART-1 (MELAN A) is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
<b>Applications</b>	IHC-F/IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	M2-9E3
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Concentrate
<b>Dilution Range</b>	IHC-F IHC-P 1:100-1:200
<b>Formulation</b>	Bioreactor Concentrate with 0.05% Azide
<b>Isotype</b>	IgG2bk
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Recombinant hMART-1 protein
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	