

**Anti-STAT6 antibody [D1] (STJ16101849)**

STJ16101849

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-STAT6 is suitable for use in Immunohistochemistry research applications.
<b>Applications</b>	IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	D1
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Concentrate
<b>Dilution Range</b>	IHC-P 1:100-1:200
<b>Formulation</b>	PBS Buffer, with less than 0.1% sodium azide and 0.1% gelatin
<b>Isotype</b>	IgG2b
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">6778</a>
<b>Gene Symbol</b>	<a href="#">STAT6</a>
<b>Uniprot ID</b>	<a href="#">STAT6_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	