

Anti-OCT-3/4 antibody [C-10] (STJ16101853)

STJ16101853

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-OCT-3/4 is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-F/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	C-10
Concentration	
Conjugation	Unconjugated
Purification	Purified
Dilution Range	IHC-F IHC-P 1:50-1:100
Formulation	Purified antibody in PBS with 0.2 % BSA and 15mM sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Amino acids1-134 of Oct-3/4 of human origin
Immunogen Region	
Specificity	
Immunogen Sequence	