

Anti-MDM2 antibody [SMP14] (STJ16101854)

STJ16101854

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MDM2 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SMP14
Concentration	
Conjugation	Unconjugated
Purification	Purified
Dilution Range	IHC-F
	IHC-P 1:50-1:100
	WB
Formulation	Affinity purified antibody in aqueous buffered solution containing $\leq 0.09\%$ sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4193
Gene Symbol	MDM2
Uniprot ID	MDM2_HUMAN
Immunogen	Human MDM2 aa. 154-167
Immunogen Region	
Specificity	
Immunogen Sequence	