

Anti-SOX10 antibody [SOX10/1074] (STJ16101855)

STJ16101855

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-SOX10 is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-F/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SOX10/1074
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-F
	IHC-P 1:200-1:400
Formulation	Bioreactor Concentrate with 0.05% Azide
Isotype	IgG2bk
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6663
Gene Symbol	SOX10
Uniprot ID	SOX10_HUMAN
Immunogen	Recombinant fragment (155 Amino residues between aa 100-300) of human SOX10 protein
Immunogen Region	
Specificity	
Immunogen Sequence	