

**Anti-CD43 antibody [DF-1] (STJ16101856)**

STJ16101856

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD43 is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC-P/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	DF-1
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Concentrate
<b>Dilution Range</b>	IHC-P 1:100-1:200 WB 1:100-1:200
<b>Formulation</b>	Bioreactor Concentrate with 0.05% Azide
<b>Isotype</b>	IgG2k
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	KG1 cells (myoblast cell line)
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	