

Anti-E-CADHERIN antibody [MX020] (STJ16101857)

STJ16101857

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-E-CADHERIN is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MX020
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:100-1:200
Formulation	Tissue culture supernatant containing 15mM sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence