

Anti-CD31 antibody [MX032] (STJ16101858)
STJ16101858**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD31 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MX032
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:200-1:400
Formulation	Tissue culture supernatant containing 15mM sodium azide.
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence