

Anti-PAX8 antibody [PAX8/1492] (STJ16101859)

STJ16101859

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-PAX8 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PAX8/1492
Concentration	
Conjugation	Unconjugated
Purification	Purified
Dilution Range	IHC-P 1:200
	WB
Formulation	Purified antibody in PBS with 0, 05% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7849
Gene Symbol	PAX8
Uniprot ID	PAX8_HUMAN
Immunogen	Recombinant human PAX8 protein
Immunogen Region	
Specificity	
Immunogen Sequence	