

Anti-Caldesmon (HMW) antibody [H-CALD] (STJ16101860)

STJ16101860

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Caldesmon (HMW) is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H-CALD
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:400-1:800
Formulation	Bioreactor Concentrate with 0.05% Azide
Isotype	IgGk
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Crude human uterus extract
Immunogen Region	
Specificity	
Immunogen Sequence	