

**Anti-SMAD4 (DPC4) antibody [B-8] (STJ16101862)**

STJ16101862

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-SMAD4 (DPC4) is suitable for use in Immunohistochemistry research applications.
<b>Applications</b>	IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	B-8
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purified
<b>Dilution Range</b>	IHC-P 1:200-1:400
<b>Formulation</b>	Purified antibody fraction from mouse anti-serum with 0.2% BSA and 15mM sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Amino acid 1-552 representing full length Smad4 of human origin
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	