

Anti-CD5 antibody [C5/473] (STJ16101863)

STJ16101863

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD5 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	C5/473
Concentration	
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:200-1:400
Formulation	Bioreactor Concentrate with 0.05% Azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	921
Gene Symbol	CD5
Uniprot ID	CD5_HUMAN
Immunogen	Human CD5 recombinant protein
Immunogen Region	
Specificity	
Immunogen Sequence	