

Anti-hCG (Alpha 3 Epitope) antibody [INN-hCG-5] (STJ16101896)

STJ16101896

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-hCG (Alpha 3 Epitope) is suitable for use in Immunohistochemistry and Immunohistochemistry research applications.
Applications	IHC-F/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	INN-hCG-5
Concentration	4 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	hCG
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	