

Anti-Aggrecan antibody [70] (STJ16101898)

STJ16101898

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Aggrecan is suitable for use in Immunohistochemistry, ELISA, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F/ELISA/IP/WB
Host/Source	Mouse
Reactivity	Human/Bovine

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	70
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Aggrecan
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	