

Anti-Beta Endorphin antibody [B31.15] (STJ16101900)

STJ16101900

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta Endorphin is suitable for use in Immunohistochemistry, Immunofluorescence and Immunohistochemistry research applications.
Applications	IHC-F/IF/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B31.15
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human beta endorphin (full length native protein)
Immunogen Region	
Specificity	
Immunogen Sequence	