

**Anti-CA125 antibody [X325] (STJ16101901)**

STJ16101901

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CA125 is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	ELISA/IHC-P/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	X325
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purified
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with <0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	CA125 antigen purified from human ovarian carcinoma
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	