

Anti-CD10 antibody [SN5c] (STJ16101902)

STJ16101902

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD10 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F/FC/IP/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SN5c
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Partially purified cell membrane antigens from fresh leukemia cells
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	