

Anti-CD161 antibody [B199.2] (STJ16101907)

STJ16101907

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD161 is suitable for use in Immunohistochemistry, Flow Cytometry and Immunoprecipitation research applications.
Applications	IHC-F/FC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B199.2
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified human NK cells cultured in IL-2
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	