

**Anti-CD29 antibody [4B7R] (STJ16101910)**

STJ16101910

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD29 is suitable for use in Flow Cytometry, Immunofluorescence, Immunoprecipitation and Immunohistochemistry research applications.
<b>Applications</b>	FC/IF/IP/IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	4B7R
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purified
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with <0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Ocular melanoma cell line V+B2
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	